


Skin Cancer Classification


User Guide – How to Use System

Disclaimer Page

The user will be prompted with a disclaimer page. To proceed with the application, the user must acknowledge the message and press the button.




Classification Page



1. To access records, click the “RECORDS” button
2. The user must upload an image prior to classification. Image types that are supported are the following:
 - a. PNG
 - b. JPEG
 - c. HEIC
 - d. WEBP
 - e. TIFF
3. After uploading an image, the user run the classification with the “RUN” button
4. The “Results” button will display the classification results
5. The “Info” button will display a reference image and additional information
6. The “Links” button will display further reference links
7. The top left button will redirect to the help page


Records Page



1. The “ANALYZE” button will redirect to the classification page
2. The “Records” button display all records
3. The “Stats” button will display a bar chat count of all the records
4. Image will display after clicking on a record
5. Right clicking a record will allow the user to delete a single or all records

Stats Page

1. The user can go back to records with the bottom button
2. Hovering over a bar will display further information



1

Help Page

The user can return to the classification page with the top right button

